Description

Multilink® Automix is a self-curing luting composite with light-curing option for the adhesive luting of indirect restorations made of metal, metal-ceramic, all-ceramic and composite. Multilink Automix is exclusively applied together with the self-etching and self-curing Multilink Primer. Multilink Primer is used to condition the dental hard tissue and any existing composite build-up fillings.

The special filler composition gives Multilink Automix a very high radiopacity. Monobond Plus is recommended as a coupling agent to achieve a strong bond to precious and non-precious alloys, as well as to all-ceramics made of zirconium and aluminium oxide and silicate ceramics.

Shades

Multilink Automix is available in four shades with different degrees of translucency:

- transparent (high translucency)
- yellow (high translucency)
- opaque (low translucency)
- white (medium translucency)

Working time

The working and setting times depend on the ambient temperature. Once Multilink Automix has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C ± 1 °C (73 °F ± 2 °F)	Intraorally (in combination with Multilink Primer) approx. 37 °C ± 1 °C (99 °F ± 2 °F)
Working time	approx. 3 min	approx. 2 min
Curing time (including working time)	approx. 8 min	approx. 5 min

Mixing ratio

Multilink Automix is always dispensed from the automix syringe in the optimum ratio. Multilink Primer A and Multilink Primer B are mixed in a 1:1 ratio (e.g. 1 drop Primer A + 1 drop Primer B).

Composition*Multilink Automix*

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers include barium glass, ytterbium trifluoride and spheroid mixed oxide. The particle size is 0.25–3.0 µm. The mean particle size measures 0.9 µm. The total volume of inorganic fillers is approximately 40 %.

Multilink Primer A and B

Multilink Primer A is an aqueous solution of initiators. Multilink Primer B contains HEMA, phosphonic acid and methacrylate monomers.

Multilink® Automix

Multilink® Automix

EN	Instructions for Use Self-curing luting composite with light-curing option Self-curing and self-etching primer	ES	Instrucciones de uso Composite de cementación universal autopolimerizables con opción de fotopolimerizable Primer autopolimerizable y de autograbado
DE	Gebrauchsinformation Selbsthärtendes zahnärztliches Befestigungskomposit mit optionaler Lichthärtung Selbsthärtender und selbstätzender Primer	PT	Instruções de Uso Compósito de cimentação autopolimerizável, com fotoativação opcional Primer autopolimerizável e autocondicionante
FR	Mode d'emploi Composite autopolymérisant avec option photopolymérisation Primer autopolymérisant et automordançant	SV	Bruksanvisning Självhårdande kompositcement för dentalt bruk, som kan ljushärdas Självhårdande och självetsande primer
IT	Istruzioni d'uso Cemento composito dentale autoindurente con fotoindurimento opzionale Primer autoindurente ed automordenzante	EL	Οδηγίες Χρήσεως Αυτοπολυμεριζόμενι ρητλώδης κοιλία με δυνατότητα φωτοπολυμερισμού Αυτοπολυμεριζόμενο και αυτοαδροποιούμενο primer
FI	Käyttöohjeet Kemiallisokvetiteinen yhdistelmämuoviselementti, jossa on valokovetusmahdollisuus. Kemiallisokvetiteinen ja itse-etsaava primeri	TR	Kullanma Talimatları Kendiliğinden sertleşen, ışıkla sertleşme opsiyonu olan yapıştırma kompoziti Kendiliğinden sertleşen ve kendiliğinden asitleyen (self-etching) primer
NO	Bruksanvisning Selvherdende dental sementeringskomposit med valgfri lysberding Selvherdende og selvetsende primer	RU	Инструкция Самотвердеющий стоматологический фиксирующий композиционный цемент с опциональной возможностью фотополимеризации Самотвердеющий самопротравливающий праймер
NL	Gebruiksaanwijzing Zelfuithardend tandheelkundig bevestigingscomposit met optie voor lichtuitharding Zelfuithardende en zelfetsende primer	PL	Instrukcja stosowania Samopolimeryzujący, kompozytowy cement adhezyjny o podwojnym mechanizmie polimeryzacji Samowytwarzający chemoutwardzalny Primer

Indication

Multilink Automix and Multilink Primer are used for the permanent cementation of indirect restorations where a strong bond is desired:

- Inlays, onlays, crowns, bridges and root posts made of
 - metal and metal-ceramics
 - all-ceramics, in particular opaque zirconium oxide ceramics
 - composites and fibre-reinforced composites

Contraindication

The use of Multilink Automix and Multilink Primer is contraindicated

- if a dry working field cannot be established or the stipulated working technique cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Multilink Automix and Multilink Primer A/B.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to single components have been reported.

Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided. Alkaline jet media impair the effect of Multilink Primer.

Hydrogen peroxide and other disinfectants with an oxidative effect may interact with the initiator system, which in turn can impair the curing process. Therefore, the automix syringe must not be disinfected with oxidative agents. The syringe can be disinfected, e.g. by wiping with medical alcohol.

Application

1. Removal of the temporary restoration and cleaning of the cavity

Remove any residue of the temporary luting cement from the cavity or preparation with a polishing brush and an oil and fluoride-free cleaning paste (e.g. fluoride-free Proxyt). Rinse with water spray. Subsequently, dry with air that is free of water and oil. Avoid overdrying.

Note: Cleaning with alcohol could lead to dehydration of dentin.

2. Try-in of the restoration and isolation

Seat the restoration using the desired Multilink Automix Try In paste to check the shade effect, accuracy of fit and occlusion of the restoration.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture. If necessary, make adjustments with fine diamonds at medium speed and with slight pressure using copious amounts of water. Polish ground surfaces. Ensure appropriate relative or complete isolation using auxiliary devices such as OptraGate® or OptraDam®.

If the dental hard tissue is contaminated with blood or saliva during the try-in procedure, it must be cleaned again as described in Item 1.

3. Pre-treatment of the restoration

3.1 If a restoration that has been pre-treated in the dental lab comes in contact with blood or saliva during the try-in procedure, the contaminated bonding surface of the restoration must be cleaned in the following way:

- After try-in, thoroughly rinse the restoration with water spray and dry it with oil-free air.
- Completely coat the bonding surface of the restoration with Ivoclean using a microbrush or brush.
- Allow Ivoclean to react for 20 seconds. Thoroughly rinse the restoration with water spray and dry it with oil-free air.
- Use Monobond Plus on restorations made of lithium disilicate glass-ceramics or zirconium oxide ceramics (e.g. IPS e.max® Press/CAD, IPS e.max ZirCAD). Apply Monobond Plus to the cleaned surfaces using a brush or microbrush. Allow Monobond Plus to react for 60 seconds and subsequently disperse it with a strong stream of air.

3.2 Restorations that have not been pre-treated by the dental lab should be rinsed with water spray and dried after the try-in procedure.

3.2.1 Thereafter, the restoration is conditioned as follows:

- a) Glass-ceramic restorations (e.g. IPS Empress®)
 - Etch the restoration with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 seconds or according to the instructions of the manufacturer of the restorative material.
 - Thoroughly rinse the restoration with water spray and dry it with oil-free air.
- b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max Press, IPS e.max CAD)
 - Etch the restoration with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 seconds or according to the instructions of the manufacturer of the restorative material.
 - Thoroughly rinse the restoration with water spray and dry it with oil-free air.
- c) Zirconium oxide (e.g. IPS e.max ZirCAD) and aluminium oxide ceramic restorations
 - Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
 - If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
 - Thoroughly rinse the restoration with water spray and dry it with oil-free air.
 - **IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.
- d) Metal or metal-supported restorations
 - Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material) until a uniformly mat surface is achieved.
 - If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
 - Thoroughly rinse the restoration with water spray and dry it with oil-free air.
 - **IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

3.2.2 Next, apply Monobond Plus with a brush or microbrush to the pre-treated surfaces. Allow the material to react for 60 seconds. Disperse it with a strong stream of air.

3.3 In general, restorations made of composite or fibre-reinforced composite must be conditioned according to the directions of the manufacturer of the respective restorative material.

4. Mixing of Multilink Primer A and Multilink Primer B

The two primer liquids Multilink Primer A and B are mixed in a 1:1 ratio (e.g. 1 drop Primer A and 1 drop Primer B). The mixed Primer A/B is solely self-curing and does not need to be protected against light. However, it must be applied within 10 minutes.

5. Application of the mixed Multilink Primer A/B to the enamel, the dentin and the build-up.

Apply the mixed Multilink Primer A/B onto the entire bonding surface using a microbrush, starting with the enamel surface. Subsequently, scrub it into to the surface for another 30 seconds.

Disperse excess with blown air until the mobile liquid film is no longer visible. As the Primer is solely self-curing, no light-curing is necessary!

Note: Previous conditioning of enamel with phosphoric acid is required for Maryland bridges.

6. Application of Multilink Automix to the restoration

For each application, place a new automix tip on the syringe. Dispense Multilink Automix from the automix syringe and apply the desired amount onto the restoration. As the luting material will cure in the used mixing tip, it serves as a seal for the remaining contents of the syringe until needed again (replace with a new tip before the next use).

Note: Multilink Automix should be used quickly after it has been dispensed and the restoration should be placed rapidly! Multilink Automix is applied directly to the inner surface of the restoration. The application of Multilink Automix directly to the preparation or the cavity, which has been pre-treated with Multilink Primer, is not indicated as this would lead to considerable acceleration of the curing process and thus may compromise the fit of the restoration.

7. Placement of the restoration and removal of excess cement

a) solely self-curing

Seat the restoration in place and secure it. Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or scaler.

Remove excess rapidly in difficult-to-reach areas (proximal areas, gingival margins, pontics). Due to the reaction between Multilink Automix and Multilink Primer A/B, a high bond strength and a high degree of cure is achieved within a few minutes of seating the restoration.

b) self-curing with additional light-curing (quarter technique, indicated for cases with up to 2 bridge abutments = 3 to 4-unit bridges)

Light-curing of all excess cement and subsequent removal: Seat the restoration in place and secure it. Light-cure excess cement in quarter segments (mesio-oral, disto-oral, mesio-buccal, disto-buccal) using a curing light at a distance of max. 10 mm. In the process, the following parameters must be observed:

Light intensity	Exposure time per quarter segment	
approx. 650 mW/cm ²	3 s	e.g. Bluephase in LOW POWER mode
approx. 1,000 mW/cm ²	1–2 s	e.g. Bluephase Style or Bluephase in HIGH POWER mode

ivoclar
vivadent[®]
clinical

CE 0123

ivoclar
vivadent[®]
clinical

Rx ONLY

Date information prepared:

2018-11-16/Rev. 3

645942/WE3

Manufacturer

Ivoclar Vivadent AG

Bendererstrasse 2

9494 Schaan/Liechtenstein

www.ivoclarvivadent.com

Thereafter, excess cement is easy to remove with a scaler. Excess must be removed rapidly in difficult-to-reach areas (proximal areas, gingival margins, pontics). Subsequently, light-cure all margins again for 20 seconds (approx. 1,000 mW/cm², e.g. Bluephase in HIGH POWER mode or Bluephase Style).

c) self-curing with additional light-curing (quarter technique, indicated for cases with up to 6 bridge abutments = circular bridges)

Light-curing of all excess cement and subsequent removal: Seat the restoration in place and secure it. Successively light-cure excess cement on all the abutments in quarter segments (mesio-oral, disto-oral, mesio-buccal, disto-buccal) using a curing light at a distance of max. 10 mm. In the process, the following parameters must be observed:

Light intensity	Exposure time per quarter segment	
approx. 650 mW/cm ²	1 s	e.g. Bluephase in LOW POWER mode

Thereafter, excess cement is easy to remove with a scaler. Excess must be removed rapidly in difficult-to-reach areas (proximal areas, gingival margins, pontics). Subsequently, light-cure all the margins again for 20 seconds (approx. 1,000 mW/cm², e.g. Bluephase in HIGH POWER mode or Bluephase Style).

Note

Like all composites, Multilink Automix is subject to oxygen inhibition: In other words, the surface layer (approx. 100 µm) does not polymerize during the curing process, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

8. Finishing of the restoration

- Finish the cement joints with preparation diamonds if necessary.
- Polish restoration margins/cement joints with suitable polishers (e.g. Astropol®) as well as finishing and polishing strips.
- Check occlusion and function and make adjustments if necessary.

Special instructions for the cementation of endodontic posts

1. For the cementation of endodontic posts, carefully clean the root canal to remove any residue of root canal fillers. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite). Ideally, the mixed Multilink Primer A/B should be applied to the root canal and the occlusal surface of the preparation for about 15 seconds each using the thin (violet) microbrushes. Remove excess from the root canal using paper points.
2. Coat the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed Multilink Automix cement.
IMPORTANT! Do not apply Multilink Automix cement with the lentulo into a root canal that has been wetted with Primer A/B. Otherwise, premature curing might occur, which could possibly prevent the post from being properly positioned.
3. Place the endodontic post. Excess cement will be displaced.
4. Distribute excess cement across the occlusal preparation surface for total coverage. Then, light-cure Multilink Automix for 20 seconds.
5. Core buildup material (e.g. MultiCore®) can then be directly applied onto Multilink Automix to serve as a bonding agent and cured according to the instructions of the manufacturer (light-curing).

Warning

Unpolymerized Multilink Automix pastes and the Multilink Primer liquids may cause slight irritation. Avoid contact with the skin, mucous membrane and eyes. If the material comes into contact with the eyes, immediately rinse with plenty of water and seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Multilink Automix and Multilink Primer must not be used after the expiry date.
- Storage temperature 2 – 28 °C / 36 – 82 °F.
- Multilink Primer bottles must be tightly closed with the correct cap after use.
- In order to ensure a tight seal of Multilink Automix syringes, the mixing tip is left on the syringe after use.
- Expiry date: See information on bottles, syringes and packaging.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.